

# Measurement of pH and Titratable Acidity

D. Garner, C.H. Crisosto, P. Wiley, and G.M. Crisosto

## I. Materials

**A. Required:** pH meter or phenolphthalein, burette, burette clamp and stand, gram scale, graduated cylinder, beakers, 0.1N NaOH solution

**B. Optional:** magnetic stirrer & stir bar, automatic titrator

## II. Procedure

**A. Obtain at least 50 mls of clear juice by one of the following methods:**

1. Cut fruit, press with a hand press, and filter through cheesecloth, or
2. Cut fruit into a blender, homogenize, centrifuge slurry, and pour off clear liquid for analysis.

\*\* Sugar levels often vary within the fruit, being higher at the stem-end and lower at the calyx-end. For this reason, it is important to use longitudinal slices of fruit (from end to end) when sampling.

**B. Make sure samples are at room temperature before taking measurements.**

**C. Measure the pH of the samples with a pH meter and record the value.**

**D. For each sample, weigh out 6 grams of juice into a 100 ml beaker.**

**E. To each sample, add 50 mls of water.**

**F. Titrate each sample with 0.1 N NaOH to an end point of 8.2 (measured with the pH meter or phenolphthalein indicator) and record the milliliters (mls) of NaOH used.**

**G. Calculate the titratable acidity using the following formula:**

$$\% \text{ acid} = \frac{[\text{mls NaOH used}] \times [0.1 \text{ N NaOH}] \times [\text{milliequivalent factor}] \times [100]}{\text{grams of sample}}$$

| Commodity                      | Predominant Acid | Milliequivalent Factor |
|--------------------------------|------------------|------------------------|
| Stone fruit, apples, kiwifruit | Malic Acid       | 0.067                  |
| Citrus                         | Citric Acid      | 0.064                  |
| Grapes                         | Tartaric Acid    | 0.075                  |